


Take Control. Stay Informed.

A simple glance to track trends in your diabetes care and effortlessly share your Omnipod DASH[®] System data to expand your care team.

Automatic Uploads.

Featured on the Omnipod DASH[®] System.



Connect PDM to Wi-Fi

Patient data
automatically uploaded

HCPs may view and interpret
data from anywhere

Username (Email)

Password

ProConnect Code



PLEASE NOTE: if your doctor has uploaded your current Omnipod® PDM to Glooko before, you already have a Glooko account. If you can't remember your login credentials, call Glooko Support at **1-800-206-6601**. They'll provide it to you and ensure that your account is connected to your doctor's.

Create your Glooko Account:

- 1 Navigate to **my.glooko.com**
- 2 Log in or click **Sign Up for Glooko Here** and fill out your information
 - Enter a ProConnect Code to automatically share your Omnipod DASH® System data with your healthcare provider
- 3 Select **I am an Omnipod® System user**



Enable automatic upload:

- 1 Navigate to settings and select **Connect to DASH® Cloud**
- 2 Log in to **PodderCentral™**
- 3 Follow the on-screen prompts

You are all set!



DON'T FORGET TO...

- Add the ProConnect code so your physician can view your data.
- Make sure that your Wi-Fi on the PDM is turned on.

If you are experiencing any issues with Glooko, please contact Glooko Support:
Call: 1-800-206-6601 **Text:** 650-720-5310
Email: support@glooko.com

For more information on the Omnipod DASH® System, contact the Insulet Customer Care Team:
Call: 1-800.591.3455 **Visit:** omnipod.com